

## **Agilia** ProNeo

Nutrition parameters

#### Nutrition flow rate range

0.1 - 200 mL/h, depending on the syringe capacity.

#### Flow rate accuracy +/- 1% on mechanism; +/- 2 % with syringes

## Syringe capacities

Types of syringe

## Nutrition mode

- Simple rate: Nutrition with a programmed rate
   Volume/Time: Nutrition of a programmed volume (0.1-99.9mL) over a programmed period of time (1min 96h) Volume Limit: Nutrition with a limitation on the volume to be delivered (0.1-999mL)

#### Volume delivered

Rate: depending on the syringe capacity (250-1200 mL/h).

Bolus rate depending on the syringe capacity (50-500 mL/h).

#### End administration (V/T & VL)

KTO (Keep Tube Open): adjustable from 0.1 to 5 mL/h, continuous nutrition or stop.

Programmable from 1 minute to 24 hours, increments minute by minute.

**Data log event** 1500 data log events in real time.

The night mode decreases the brightness of the screen and the green lights. The key beep can optionally be turned off. The night mode can be programmed manually or automatically in a variable

#### Pressure management

2 pressure controlling modes available: variable or 3 pre-set levels - Range from 50 to 900 mmHg. (25 mmHg increment from 50 to 250 mmHg / 50 mmHg increment from 250 to 900 mmHg). Can be enabled / disabled and adjusted

The Dynamic Pressure System - DPS - warns of pressure variations. A risk of obstruction or a possible leak in the nutrition line can thus be anticipated.

#### Anti-bolus system

Reduces significantly bolus after occlusion release. < 0.35 mL max for a 60 mL syringe.</p>

#### Alarms / Pre-alarms / Security

Green light for nutrition in progress, orange light for low and medium priority alarms, red for high priority alarms. All alarms are expressed by means of light indicators, written words, pictograms and sound beeps

#### Syringes installation control

Syringe barrel clasp check, plunger head detection, anti-siphon system check, flange detection.

#### Nutrition control

Occlusion pressure pre-alarm, occlusion pressure alarm, end of nutrition pre-alarm, end of nutrition alarm, volume limit pre-alarm, volume limit alarm, keypad manual locking or keypad autolock, hard and soft flow rate limits, start nutrition at pause end.

#### Device control

Disengaged driving mechanism alarm, low battery pre-alarm, discharged battery alarm, battery capacity display in hours and minutes, unconfirmed programming, technical malfunction alarm (auto-test, rotation), drive system advance check, watchdog check, communication connection failure, plug-head disengagement, autolock/lock code (on Keypad).

#### Maintenance

Preventive maintenance warning.

#### Technical specifications

**Manual pusher**Protection for the ongoing nutrition thanks to "Push-Guard".

Blue graphic LCD monochrome, size 70 x 35 mm (256 x 128 pixels).

## Swinglock clamp

Versatile clamp that allows the fixation on a rail or on a pole (Pole: diameter 15 to 40 mm max. / Rail: 25 to 35 x 10 mm).

Stackability
Up to 3 devices self-stackable on a pole.

## Dimensions (h/w/d) / weight

**Battery**Characteristics: 7.2 V 2.2 Ah - Li-ion Smart battery, remaining battery life and battery charge level available on the display.

Battery Life (when fully charged): > 11 h at 5 mL/h.

Battery recharge: - Pump OFF: < 6 h - Pump ON: < 20 h

## Waterproofness

#### Power supply

100 V - 240 V ~ / 50 / 60 Hz with functional earth

### Compliance

## Electromagnetic compatibility EMC

## Safety & essential performance of pump

CE 0123 marking in compliance with the Council Directive 93/42/EEC

## **Electrical Compliance**

Medical Device Directive

Protection against leakage current: Defibrillation-proof type CF Protection against electric shocks: class II in accordance with IEC 60601-1

## Alarm system

Usability Engineering IEC 60601-1-6 and IEC 62366

Agilia and ProNeo are registered trademarks by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information.











Enteral Nutrition Syringe Pump







# **Agilia** ProNeo

# Enteral Nutrition Pump for Neonates

Agilia ProNeo is a syringe pump dedicated to Enteral Nutrition delivery for neonates, prematures babies and pediatric patients.

The pump and its disposables are designed to improve safety in neonatal and pediatric care units.



# With Agilia ProNeo, little patients are protected to grow

## Safety

Agilia ProNeo prevents possible confusion between intravenous and enteral delivery

- Purple color and labelling for specific enteral use
- Automatic detection of syringe capacity
- Pushguard and keypad lock for a secured usage





## Accuracy

Agilia ProNeo ensures respect of prescribed nutrition in neonatal intensive care unit to achieve a better outcome for neonates.

- Different delivery modes to cover clinical practice needs: Flowrate only / Volume over Time/ Volume limit
- Flowrate accuracy\* +/- 3 % for a precise milk delivery







 $\ast$  +/- 1 % on mechanism, +/- 2 % with syringes

## Convenience

Agilia ProNeo has been designed to be convenient on what matters for newborn nutrition

- Night mode for baby well-being
- Customisable pressure management for milk delivery
- Easy navigation through the menu



# ProNeo Product Line

The ProNeo product line for pump offers high flexibility to the users and reliability for the patients.



## ProNeo Feeding Tubes

Polyurethane material soft and transparent, smooth surface to ease the passage

## ProNeo Extension Sets

Polyurethane material soft and transparent

Tube clamp: prevent leakage during preparation/storage



Dedicated to enteral nutrition delivery in Neonatal and Paediatric Care Unit



Fresenius Freka Connect syringes ENFit syringes from 5mL to 60mL

