

**FRESENIUS
KABI**

caring for life



Smart Holder COM

Instructions For Use

EN	Instructions For Use	3
BG	Инструкции за употреба	22
CN	使用说明书	43
CZ	Návod k obsluze	63
DE	Gebrauchsanweisung	83
DK	Brugervejledning	104
ES	Manual de usuario	124
FI	Käyttöohjeet	144
FR	Notice d'utilisation	163
GR	Οδηγίες χρήσης	184
HR	Upute za uporabu	204
HU	Használati útmutató	225
IT	Istruzioni per l'uso	245
JP	取扱説明書	266
NL	Gebruiksaanwijzing	285
NO	Bruksanvisning	306
PL	Instrukcja obsługi	326
PT	Manual de Instruções	347
RO	Instrucțiuni de utilizare	367
RU	Руководство по эксплуатации	388
SE	Bruskanvisning	410
SK	Návod na obsluhu	430

Symbol descriptions

Symbol	Description	Symbol	Description
	Medical Device		Unique Device identifier
	Refer to the Instructions For Use		Index of protection against solid foreign objects (> 2.5 mm) and dripping liquids
	Manufacturer		Class II equipment
	Defibrillation-proof type CF applied part		Serial number
	Catalogue number		Input terminal – connector
	Output terminal – connector		Alternating current (AC)
	Direct current (DC)		CE Marking
	Part included in a recycling process		Humidity limitation
	Fragile, handle with care		Atmospheric pressure limitation
	This way up		General symbol for recyclable material
	Keep dry		Medical electrical system weight(kg)
	Temperature limit		Name and address of the manufacturing facility
	UDI (01) 04052682037238 (21) 12345678 (11) 200316 (240) CS1000429		(01) Product identifier GTIN (21) Product Serial Number (11) Date of Manufacture (240) Product Reference

Caution symbol: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

Information symbol: Recommendations to be followed.

INFORMATION

- Please refer to the Use environment section for additional information on temperature, pressure and humidity limitations, and altitude.
- These Instructions For Use (IFU) are also available online. For further information, check IFU for Amika, Amika+ or Ambix nova on <https://key2.fresenius-kabi.com>.

EN Table of contents

1	INTRODUCTION	6
1.1	Scope	6
1.2	Intended use.....	6
1.3	Use environment	6
2	DESCRIPTION	7
2.1	System definition	7
2.2	Packaging content	7
2.3	Holder description.....	7
3	HOLDER INSTALLATION AND REMOVAL	8
3.1	Installation.....	8
3.2	Removal	9
4	NURSE CALL CONNECTION AND REMOVAL	9
4.1	Specifications and use.....	10
4.2	Connection to the holder	10
5	NURSE CALL QUICK CHECK PROTOCOL	11
6	RS232 SERIAL COMMUNICATION	12
6.1	Communication with PDMS	12
6.2	USB Maintenance cable connection and removal.....	12
7	TROUBLESHOOTING	13
8	TECHNICAL INFORMATION	14
8.1	Operation mode.....	14
8.2	Power supply specifications	14
8.3	Nurse Call Specifications.....	14
8.4	RS232 serial communication connector specifications.....	14
8.5	Dimensions – Weight	15
8.6	Compliance with standards	15
9	CLEANING AND DISINFECTING	15
9.1	Prohibited cleaning or disinfection agents	15
9.2	Precautions	15
9.3	Recommended cleaning and disinfection agents	16
9.4	Cleaning and disinfection guidelines and protocol	16
10	TRANSPORT, STORAGE AND RECYCLING CONDITIONS	17
10.1	Storage and transport conditions.....	17
10.2	Storage.....	17
10.3	Recycling and Disposal	18
11	GUIDANCE AND MANUFACTURER'S DECLARATION OF EMC	18

12 SERVICES	18
12.1 Warranty	18
12.2 Limited warranty	18
12.3 Maintenance requirements.....	19
12.4 Service policy and rules.....	19
12.5 Notification of serious incident	19
13 ORDERING INFORMATION	19
14 GLOSSARY OF TERMS	20

1 Introduction

This document is written to the attention of:

End users of Amika / Amika+ / Ambix nova with Smart Holder COM.

1.1 Scope

These Instructions for Use (IFU) are applicable to the Smart Holder COM referred to as holder. The holder is compatible with Amika enteral feeding pump, Amika+ enteral feeding pump and Ambix nova ambulatory infusion pump, which are referred to as pump in this IFU.

CAUTION

The user must follow the instructions specified in this IFU. Failure to observe these instructions may result in damage to the equipment, injury to patients or injury to users. Specific texts are highlighted using the symbols described in *Symbol descriptions*.

1.2 Intended use

The holder can be positioned on authorized support such as a pole, rail, bed, wheelchair, table and table stand. It allows the pump to be attached to it to ensure the pump is securely positioned.

It is intended to supply the operating power and to charge the battery for the pump installed on it when it is connected to the mains power supply.

It is intended to connect a pump to an external Nurse Call system to transmit a pump alarm state.

It is intended to connect a pump to a PC for service activities via Amika and Ambix nova Partner software.

It is intended to communicate between a pump and a third party system PDMS (Patient Data Management System).

1.2.1 Intended user population

The holder must only be used by trained healthcare professionals, patients, and any other relevant users. For training guide, contact your Fresenius Kabi sales representative.

1.3 Use environment

The holder is intended to be used inside the hospital, nursing station, home, biomedical office and any other areas involved in technical service.

CAUTION

- The holder should be used under specified operational, storage and transport conditions listed below.
- DO NOT USE : in ambulance, helicopter and aircraft
in Magnetic Resonance Imaging (MRI) environments

	Operating	Storage and Transport
Temperature Range	10°C to 40°C	-20°C to +45°C
Pressure Range	700 hPa to 1060 hPa	500 hPa to 1060 hPa
Humidity Range	30% to 85%, no condensation	10% to 90%, no condensation
Altitude	Less than 3000m	-

2 Description

2.1 System definition

Place Amika / Amika+ / Ambix nova pump on the holder to form an integrated Amika / Amika+ / Ambix nova system.

2.2 Packaging content


The holder is delivered with the following:

- Smart Holder COM
- Power cable
- Nurse Call cable
- User documents

Packaging consists of: recycled cardboard.

Symbols used on packaging are described in *Symbol descriptions*.

2.3 Holder description

Legend

- ① Clamp handle
- ② Contact pins for pump to holder connection
- ③ Slot
- ④ Grey locking lever
- ⑤ Mains supply light indicator
- ⑥ Pole clamp
- ⑦ Holder identification label
- ⑧ Power cable inlet
- ⑨ Nurse call connector
- ⑩ RS232 serial communication connector

Note: The slot can be used to fix the cables compatible with its dimension.

3 Holder installation and removal


CAUTION

Give particular attention to the risk of strangulation with cables and sets, and with the small parts that could be swallowed or inhaled.

3.1 Installation



3.1.1 Using the pole clamp

The holder can be attached universally, vertically and horizontally. Turn the pole clamp to the suitable position.

3.1.2 Positioning the holder on a rail, pole, bed or wheelchair

Ensure the holder is positioned so that the display is at the suitable height to ensure good visibility and orientation in the reading direction (the contact pins are at the bottom).


X, Y min = 10 mm
X, Y max = 35 mm

Ø min = 8 mm
Ø max = 40 mm

Fasten pole clamp firmly on the pole or rail to avoid any movement of the pump. Ensure that the pump is securely attached and firmly positioned.

3.1.3 Positioning the holder on a table


The holder can be placed on a flat and horizontal table as indicated in the figure.

Ensure the holder is positioned away from table edges to avoid being accidentally pushed off the table.

3.1.4 Positioning the pump

Slide the pump down until the grey locking lever locks the position.

3.1.5 Electrical connection and removal


CAUTION

ONLY use the power cable delivered with the device. The power cable and power socket must remain accessible at all times to allow emergency power supply disconnection.

Ensure power cable is not damaged.

To charge the pump battery or to use the pump on the mains power supply:

1. Connect the power cable to the holder.
2. Plug the power cable to the mains socket.
3. To disconnect, unplug power cable from power socket and remove power cable from holder.


The mains power supply is indicated by a green light on the holder's front panel as well as on the pump's front panel (keypad).

A beep is emitted by the pump when the power cable is disconnected.

3.1.6 Attaching and Removing the Quick Guide


A quick guide can be easily attached and removed from the pump holder as indicated on the figure of the right side.


3.2 Removal

3.2.1 Removing the pump from the pump holder

- 1 Push the grey locking lever.
- 2 Pull the pump up.

3.2.2 Removing the pump holder

4 Nurse Call connection and removal

The Nurse Call connection is functional only if:

- the pump is correctly installed on the holder,
- the holder is connected to the power supply,
- the Nurse Call cable is correctly plugged.

If the Nurse Call is not functional, pump alarm state is not transmitted.

The Nurse Call system availability and technical compliance are the responsibility of the hospital.

INFORMATION


- The Nurse Call cable is delivered with an unterminated side, which will require customization following specific requirement and installation as described below.

- All connections and disconnections must be performed by qualified and appropriately trained staff.

4.1 Specifications and use

4.1.1 Nurse Call cable specification


Wires length = 5 cm Shrinkable tube = 2 cm

The unterminated side is composed of 3 wires, connected to Nurse Call signalling system:

- R:** Stand-off contact
- T:** Working contact
- C:** Common contact

4.2 Connection to the holder

- Connect the terminated Nurse Call cable to holder Nurse Call connector.
- Pinch the cable into the slot provided for this purpose.
- Connect the Nurse Call cable unterminated side to the hospital Nurse Call system as described below.
- Check that the Nurse Call system is functional by generating alarms, Refer to *Nurse Call quick check protocol*. Ensure the pump alarm is transmitted on the connected Nurse Call system.
- To disconnect, unplug the Nurse Call cable.


4.2.1 Nurse Call system electrical schematics

Without pump alarms:

System operation	Contact	Electrical continuity
Pump is running without any alarm Pump is OFF No pump is installed on the holder Pump is booting	C-T	C-COM-NO-T

With pump alarms:

System operation	Contact	Electrical continuity
Pump is running with alarm		
Pump is running with prior information to alarm	C-R	C-COM-NC-R

5 Nurse Call quick check protocol

The following protocol provides the user with a quick check guide to allow users to confirm device behaviour according to these instructions for use and to ensure the Nurse Call is functional.

CAUTION

It is recommended that this protocol be applied regularly after installation (especially before use on a new patient, after cleaning or after moving) and shall be performed by qualified and appropriately trained staff.

Location	Ward name:	Bed name:
Smart Holder COM serial number		

Test description	Expected results	Compliance	
		Conform	Not conform
Visual inspection	Nurse Call cable, holder and power cable are not altered or damaged in any way	<input type="checkbox"/>	<input type="checkbox"/>
	Check correct customization of the cable	<input type="checkbox"/>	<input type="checkbox"/>
	Check power cable and Nurse Call cable are connected and operational	<input type="checkbox"/>	<input type="checkbox"/>
Switch pump ON and generate an alarm (eg: start the pump with no giving set installed)	A visual and audible alarm is displayed on the pump	<input type="checkbox"/>	<input type="checkbox"/>
	The LED is blinking in yellow on the pump	<input type="checkbox"/>	<input type="checkbox"/>
	The pump alarm is transmitted on connected Nurse Call system	<input type="checkbox"/>	<input type="checkbox"/>

The accessory is declared functional ONLY if all tests do CONFORM.

Note: if one or more tests do not conform, please contact your Fresenius Kabi sales representative.

<p>Operator name and signature</p> <p>Date: ___ / ___ / ___</p>

6 RS232 serial communication

INFORMATION


- RS232 serial communication connector is only used during service activities with Amika and Ambix nova Partner software and during communication with PDMS.
- All connections and disconnections must be performed by qualified and trained personnel.
- All IT devices (including computers, hubs and switches) inside the patient area (<1.5m) must comply with IEC/EN 60601-1 (leakage current).
- IT devices connected outside the patient area (>1.5m) must be at least IEC/EN 60950 compliant.
- Do not disconnect communication cables while data is being transferred.

6.1 Communication with PDMS

INFORMATION

- Before connecting the pump to a hospital information system (PDMS), please contact your IT or biomedical department.
- Fresenius Kabi is not responsible of providing serial communication cable. If need, please contact your PDMS provider for serial communication cable in accordance with their specifications (DB9, USB...).
- The serial interface will automatically operational when the pump is installed on the Smart Holder COM and running.

6.2 USB Maintenance cable connection and removal

INFORMATION

- USB Maintenance cable is only used to connect the pump with Amika and Ambix nova Partner software for service activities.
- USB Maintenance cable is not provided within the package and must be ordered separately.
- Use ONLY cables recommended by Fresenius Kabi. Refer to per pump IFU, section Ordering Information.
- Before any maintenance operations, the pump must be switched off and should not be connected to the patient.
- For configuration of IT-network and technical specifications, please refer to the Amika and Ambix nova Partner User Manual (ref. 7125).

1. Connect the terminated USB Maintenance cable on holder's RS232 serial communication connector.
2. Connect the other end of the USB Maintenance cable on the PC USB connector.
3. Check connection status in the computer.
4. To disconnect, unplug the USB Maintenance cable on holder side and PC side.

7 Troubleshooting

This section provides information to perform a first troubleshooting of the holder.

INFORMATION

For additional information about issues and troubleshooting please refer to pump Technical Manual or contact your Fresenius Kabi sales representative.

Issue descriptions	Recommended actions
Power cable cannot be connected or removed from holder	<ul style="list-style-type: none"> ● Check the power cable inlet on holder ● Check the power cable connector ● Consult technical manual ● Contact your Fresenius Kabi sales representative
Mains supply light indicator do not work	<ul style="list-style-type: none"> ● Check the power cable connection ● Check the power cable ● Consult technical manual ● Contact your Fresenius Kabi sales representative
The holder cannot be installed on the position	<ul style="list-style-type: none"> ● Check the installed position ● Check the pole clamp ● Consult technical manual ● Contact your Fresenius Kabi sales representative
Nurse Call cable cannot be connected or removed from holder	<ul style="list-style-type: none"> ● Check cable connector ● Consult pump technical manual ● Contact your Fresenius Kabi sales representative
The Nurse Call system does not replicate pump alarms	<ul style="list-style-type: none"> ● Check Nurse Call cable connection ● Check holder is connected to the mains power supply ● Check pump and holder connectors ● Consult pump technical manual ● Contact your Fresenius Kabi sales representative
PDMS or USB connection is not functional	<ul style="list-style-type: none"> ● Check the communication cable connection ● Check the communication cable ● Check holder is connected to the mains power supply ● Check pump and holder connectors ● Consult pump technical manual ● Contact your Fresenius Kabi sales representative

8 Technical information

8.1 Operation mode

The holder is a reusable device and the mode of operation is continuous.

8.2 Power supply specifications

The power cable must be connected directly to the mains power socket.

Protection against electric shocks: class II

Holder input	AC input voltage: 100-240 Vac AC input frequency: 50/60 Hz AC input current: 110mA-205mA
Holder output	9 Vdc \pm 5 % / 9W (maximum load)
Power cable length	Approx. 2 m (except plug type M is approx. 3 m)

8.3 Nurse Call Specifications

Nurse Call Connector Specifications:

Input/output	24 Vdc SELV(Safety Extra Low Voltage)/ 0.5 A maximum 24 Vac / 0.5 A maximum
Electrical insulation	1.5 kV insulation

The following table lists the minimum requirements that the connected Nurse Call system must fulfil:

Nurse Call system specification	Minimum requirement
System type	Electrical analogical system.
Compatibility	Compatible with 3 wires contacts. Relay type, dry contacts stand-off / work contact – positive / negative polarity option.
Maximum power	DC 24V – 0.5 A maximum or AC 24V – 0.5 A maximum.
Electrical insulation	1.5 kV insulation.
Compliance with standards	IEC/EN 60601-1 (leakage current).

8.4 RS232 serial communication connector specifications

Input/output	RS232 signal
Electrical insulation	1.5 kV insulation
Compliance with standards	IEC/EN 60601-1 (leakage current)

CAUTION

- This connector is for data communication and maintenance use only. Improper use of the RS232 Serial communication connector may cause the impossibility to perform maintenance or data communication.

8.5 Dimensions – Weight

	Weight	Dimensions (H × W × D)
Holder	~ 0.45 kg	Approx. 132 mm × 118 mm × 46 mm (without pole clamp)
Nurse Call cable	~ 0.05 kg	Length: approx. 2m
Packaging	~ 0.20 kg	-

8.6 Compliance with standards

	Conform to Regulation (EU) 2017/745	Protection against moisture: Holder: IP32 (drip-protected)
Safety of ElectroMedical Equipments	Conform to EN/IEC 60601-1 EN/IEC 60601-1-6 EN/IEC 60601-1-11	 Defibrillation-proof type CF applied part
EMC	Conform to EN/IEC 60601-1-2	 Class II equipment
Labeling	Conform to EN 1041 EN ISO/ISO 15223-1	

Note: The full list of applicable standards is available upon request. The device is protected against leakage current and do not disturb ECG or EEG devices.

9 Cleaning and disinfecting

CAUTION

- Do not use the holder and return the device to your biomedical department in case contamination remains visible on the device after cleaning.
- Please follow the disinfecting and cleaning best practices to limit risks of holder damage.
- No cleaning or disinfection activities should be performed when device is connected to the mains power supply.

9.1 Prohibited cleaning or disinfection agents

Do not use cleaning or disinfection agents that contain the following substances as these aggressive agents may damage the plastic parts of the device and cause the device to malfunction:

- trichloroethylene
- abrasive detergents

9.2 Precautions

Clean the holder as soon as it becomes contaminated with tube feed or drugs, and at least once a week. After cleaning, the holder should be left to dry for approximately 5 minutes before being reconnected to the mains.

The holder must be cleaned after each patient usage by a trained nurse or assistant nurse.

9.3 Recommended cleaning and disinfection agents

Didecyldimethylammonium chloride (example: Wip'Anios Excel by Anios).

Please contact the appropriate service, responsible for cleaning and disinfection products in your establishment for further details.

9.4 Cleaning and disinfection guidelines and protocol

INFORMATION

- Do not immerse holder in liquids or let liquids enter device's housing.
- Holder is resistant to recommended cleaning agents.

9.4.1 Cleaning Instructions

Prerequisites

- All cables are unplugged.
- The pump is disconnected from the holder.
- The air is at room temperature (20 °C to 25 °C).
- The operator is wearing suitable protective equipment.

Protocol

1. Place the holder on a cleaned surface or disposable underlay.
2. During cleaning, do not turn the holder over to avoid liquid leak in.
3. Use a ready-to-use wipe to remove any major grime.
4. Thoroughly wipe down all exposed surfaces (housing, pole clamp screw, connection area, cables etc.) of the holder, from top to bottom.

A minimum cleaning of 1 minute is recommended (allow to remain visibly wet for 1 minute), until all organic matter is dissolved and removed. Do not allow liquids to run, leak, or drip into the holder housing. Use cotton wool to clean the contact pins.

5. Using a fresh ready-to-use wipe, thoroughly wipe down all exposed surfaces. A minimum cleaning of 1 minute is recommended (allow to remain visibly wet for 1 minute), until all organic matter is dissolved and removed.
6. Wipe down the cables.
7. Allow the device to dry completely at room temperature.

9.4.2 Disinfection instructions

Prerequisites

- The cleaning protocol has been performed.
- All cables are unplugged.
- The pump is disconnected from the holder.
- The air is at room temperature (20 °C to 25 °C).
- The operator is wearing suitable protective equipment.

Protocol

1. Place the previously cleaned holder on a cleaned surface or disposable underlay.
2. During disinfection, do not turn the holder over to avoid liquid leak in.
3. Use a ready-to-use wipe to wipe down all exposed surfaces of the holder, making sure to cover all cracks, crevices, and hard-to-reach areas. Do not allow liquids to run, leak, or drip into the holder housing.
4. Using a fresh ready-to-use wipe, repeat steps 3. Ensure that the minimum contact

time for each step is 3 minutes for bactericide activity (surface remain visibly wet for 3 minutes). Respect the indicated contact time from the manufacturer recommendations for the required antimicrobial activity.

5. Wipe down the cables.

10 Transport, storage and recycling conditions

10.1 Storage and transport conditions

Before transport, disconnect and remove the power cable from the holder. Disconnect all other cables.

Check power cable and Nurse Call cable are operational after transport.

The holder should be used under specified storage and transport conditions (see *Use environment*) to ensure holder performance and to avoid holder malfunction.

10.2 Storage

Please make sure the holder is stored in an appropriate manner so as to avoid holder malfunctioning.

INFORMATION

- The storage area must be clean, organized and compliant with the storing conditions mentioned above.
- The holder must be handled with care during storage.

CAUTION

The holder must be cleaned and disinfected prior to storage (see *Cleaning and disinfecting*).

10.2.1 Prepare the device for storage

In order to prepare the device before storage, proceed as specified below:

1. Disconnect power cable and all other cables.
2. Remove the pump and holder from pole or rails.
3. Clean the holder (see *Cleaning and disinfecting*).
4. Handle the holder with care and store it in a compliant area.

10.2.2 Install the device after storage

INFORMATION

- We recommend that a quick check protocol is performed (refer to pump IFU, section Quick Check Protocol) when the device is installed after transport, in case of prolonged storage, or before being used on a new patient.
- To ensure Nurse Call is functional, we recommended that the *Nurse Call quick check protocol* is applied regularly after installation (especially before use on a new patient, after cleaning or after moving) and shall be performed by qualified and appropriately trained staff.

10.3 Recycling and Disposal

Devices with this label must not be disposed of with the general waste. They must be collected separately and disposed of according to local regulations.

For further information regarding waste processing regulations, contact your local Fresenius Kabi sales representative.

11 Guidance and manufacturer's declaration of EMC

CAUTION

- The holder is intended to be used with Amika, Amika+ enteral feeding pump or Ambix nova ambulatory infusion pump in the specified electromagnetic environment. The customer or the user should ensure that it is used in such an environment. For information on EMC compliance, please refer to Amika pump IFU, Amika+ pump IFU or Ambix nova pump IFU, section Guidance and manufacturer's declaration on EMC.
- The device should not be used next to other equipment. If adjacent use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
- Prolonged exposure to X-ray environments can damage the electronic components of the device. Please check pump's IFU for recommendations of safe usage in such environments.

12 Services

12.1 Warranty

12.1.1 General conditions of warranty

Fresenius Kabi guarantees that the holder is free from defects in material and workmanship during the period defined by the accepted sales conditions, except for the accessories. Accessories may have specific conditions of warranty. Please contact your Fresenius Kabi sales representative for additional information.

12.2 Limited warranty

To benefit from the materials and workmanship guarantee from our sales representative or agent authorized by Fresenius Kabi, the following conditions must be respected:

- The device must have been used according to the instructions described in this user guide and other accompanying documents;
- The device must not have been damaged when in storage or at the time of repair, and must not show signs of improper handling;
- The device must not have been altered or repaired by non-qualified personnel;
- The serial number must not have been altered, changed, or erased.

INFORMATION

- If one or more of these conditions have been violated, Fresenius Kabi will prepare an estimate for repairs covering the parts and labour required.

- When a return and/or repair of the device are required, please contact your Fresenius Kabi sales representative.

12.3 Maintenance requirements

CAUTION

- The qualified personnel must be informed if the device is dropped or if any malfunction occurs. In this case, the device must not be used. Please contact your biomedical department or Fresenius Kabi.
- When replacing components, only use Fresenius Kabi spare parts.
- No maintenance action should be performed when device is in use.

Life time of Holder: 10 years.

12.4 Service policy and rules

If the device must be sent for servicing, proceed as specified below:

1. Contact Fresenius Kabi to have packaging shipped to your facility.
2. Clean and disinfect the device to avoid of potential harm or risks to staff.
3. Pack the device in the provided package.
4. Ship the device to Fresenius Kabi.

INFORMATION

- Fresenius Kabi is not liable for loss or damage to the device during transport.
- For further information concerning the device servicing or use, please contact our sale representative.

12.5 Notification of serious incident

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority.

Information and contact information:

Fresenius Kabi AG

Else-Kröner-Str. 1

61352 Bad Homburg

GERMANY

Tel.: +49 (0) 6172 / 686-0

<http://www.fresenius-kabi.com>

13 Ordering information

CAUTION

- The holder is compatible with the range of Amika pumps, Amika+ pumps and Ambix nova pumps only.
- Use only recommended accessories from Fresenius Kabi.


For operating with these pumps or any further information, please contact our sales Department.

For further information regarding Amika / Amika+ / Ambix nova accessories, please refer to Amika pump IFU, Amika+ pump IFU and Ambix nova pump IFU, section Ordering information.

EN 14 Glossary of terms

Term	Description
A	Amper
°C	Celsius Degree
CE marking	European Conformity Mark
cm	Centimeters
ECG	Electrocardiogram
EEG	Electroencephalogram
EMC	Electromagnetic compatibility
H × W × D	Height × Width × Depth
hPa	Hecto Pascal
Hz	Hertz
IEC	International Electrotechnical Commission
kg	Kilogram
LED	Light Emitting Diode
m	Meters
mA	Milliampere
Max.	Maximum
Min.	Minimum
mm	Millimeters
MRI	Magnetic Resonance Imaging
PDMS	Patient Data Management System
V	Volt
Vac	Volt Alternating Current
Vdc	Volt Direct Current
W	Watt

Local contacts for servicing

This document may contain inaccuracies or typographical errors. Modifications may thus be made, and included in later editions. Due to the evolution of standards, and of legal texts and materials, the characteristics indicated in the text and images of this document are applicable only to the device with which it is included.

This document may not be reproduced in whole or in part without the written consent of Fresenius Kabi.

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Fresenius Kabi (Nanchang) CO., Ltd.