

Ordering information

	Size in ml	Packaging	Art.No.
Enteral Dispensing Systems, ENFit			
Freka® Connect, ENFit/ProNeo	10 ml	10 x 1	9000784XS
	10 ml	30 x 1	9000784S
	10 ml	100 x 1	9000784
	20 ml	10 x 1	9000785XS
	20 ml	30 x 1	9000785S
	20 ml	100 x 1	9000785
	60 ml	10 x 1	9000786XS
	60 ml	30 x 1	9000786S
	60 ml	50 x 1	9000786
Eccentric	60 ml	50 x 1	9000789
	100 ml	10 x 1	9000787XS
	100 ml	30 x 1	9000787
Low Dose Tip Syringes			
Freka® Connect ENFit/ProNeo LDT	1 ml	30 x 1	9000788LDT
	2.5 ml	30 x 1	9000782LDT
	5 ml	30 x 1	9000783LDT
Accessories			
Freka® Cap ENFit/ProNeo, sterile		30 x 1	7981402
Freka® Cap ENFit/ProNeo		1 x 30	7981403
ENFit Sampling Straw	15 cm	1 x 30	7752013
ProNeo Sampling Straw, sterile	10 cm	30 x 1	7752017
	15 cm	30 x 1	7752016
ProNeo Sampling Cannula, sterile	5 cm	30 x 1	7752012
ProNeo Syringe Adapter, sterile		100 x 1	7752014

The Freka® Connect Syringes are DEHP- and Latex-free.
 For more information on the ENFit syringes please contact your local Fresenius Kabi representative or visit the website: www.stayconnected.org.

© Fresenius Kabi Deutschland GmbH. Vervielfältigung und Weitergabe - auch auszugsweise - nur mit vorheriger, schriftlicher Genehmigung zulässig

739289/12 (12.22/AC)

Freka® Connect, ENFit



Freka® Connect, ENFit


Enteral Dispensing Systems, for administration of nutrition, liquids and medications

Safety advantages

- All syringes are equipped with the new safe enteral connector ENFit, according to international safety standard ISO 80369-3
- Prevention of accidental misconnections
- Tight and reliable connection of the new ENFit connector enhances the safety of the patients and caretakers

Go for safety

The ENFit syringes and the ENFit Low Dose Tip syringes reduce the risk of an unintended misconnection. Thus the ENFit and ENFit Low Dose Tip syringes offer the highest safety in the clinical practise.



ENFit =
one enteral connector
for all applications in
enteral nutrition

ENFit Low Dose Tip syringes

What are the Low Dose syringes needed for?

- The ENFit Low Dose Tip syringe is the solution for precise application and it reduces the risk of misconnection with other systems