

Accessories - Freka Connect



	Packaging	Art.-no.
Enteral Dispensing Systems, ENFit		
Freka Connect 10 ml, ENFit/ProNeo	10	9000784XS
Freka Connect 10 ml, ENFit/ProNeo	30	9000784S
Freka Connect 10 ml, ENFit/ProNeo	100	9000784
Freka Connect 20 ml, ENFit/ProNeo	10	9000785XS
Freka Connect 20 ml, ENFit/ProNeo	30	9000785S
Freka Connect 20 ml, ENFit/ProNeo	100	9000785
Freka Connect 60 ml, ENFit/ProNeo	10	9000786XS
Freka Connect 60 ml, ENFit/ProNeo	30	9000786S
Freka Connect 60 ml, ENFit/ProNeo	50	9000786
Freka Connect 100 ml, ENFit/ProNeo	10	9000787XS
Freka Connect 100 ml, ENFit/ProNeo	30	9000787

© Fresenius Kabi Deutschland GmbH. Reproduction and distribution - only with prior written authorisation. For Health Care Professionals only.



Freka® Feeding Tubes, ENFit


Freka Feeding Tubes, ENFit

- > Safe enteral connector ENFit, according to the new international safety standard ISO 80369-3
 - for a safe and tight connection
- > High kink stability
 - for a reliable flow and easy insertion
- > Integrated spiral mandrin
 - connectable to ENFit syringes

Luer-free
Latex-free
DEHP-free
PVC-free



Freka Linear Connector



Freka Tubes, ENFit are for gastric and intestinal feeding. Freka Tubes are made of polyurethane.



- Enteral connector ENFit with lean, linear connector (Fig. 1)
- Three X-ray-opaque strips embedded in the tube material for good visibility in the X-ray image
- One X-ray-radio-opaque strip embedded in the tube material with Freka Tube, FR 15, ENFit
- Tube imprinted with article number, FR-Size and tube length for easy identification (Fig. 2)
- Transparent tube material for visual inspection of administration of nutrition and liquid
- Flexible, X-ray-opaque polyurethane olive with two openings at the distal tube tip for FR 8, FR 10, FR 12 (Fig. 3)
- Tube is numbered at 5 cm steps and marked at 1 cm steps
- Product package includes: Plaster for fixing the tube

	External ø	Internal ø	Length	gastric/ intestinal	Packaging	Art.-no.
Freka® Tubes, ENFit, Linear Connector						
Freka Tube FR 8, ENFit Linear	2.8 mm	1.9 mm	80 cm	g	10 x 1	7981839
Freka Tube FR 8, ENFit Linear	2.8 mm	1.9 mm	120 cm	g, i	10 x 1	7981840
Freka Tube FR 10, ENFit Linear	3.4 mm	2.4 mm	80 cm	g	10 x 1	7981845
Freka Tube FR 10, ENFit Linear	3.4 mm	2.4 mm	120 cm	g, i	10 x 1	7981842
Freka Tube FR 12, ENFit Linear	4.0 mm	2.9 mm	120 cm	g, i	10 x 1	7981843
Freka Tube FR 15, ENFit Linear	4.5 mm	3.5 mm	120 cm	g, i	10 x 1	7981844
Freka Tube, FR 15*	4.5 mm	3.5 mm	100 cm	g	10 x 1	7980115

*for decompression, with funnel connector

gastric /
intestinal

Freka Y-Connector




Freka Tubes, ENFit are for gastric and intestinal feeding. Freka Tubes are made of polyurethane.

- Enteral connector ENFit for ENFit syringe and decompression (Fig. 1)
- Three X-ray-opaque strips embedded in the tube material for good visibility in the X-ray image
- One X-ray-opaque strip embedded in the tube material with Freka Tube, FR 15, ENFit
- Tube imprinted with article number, FR-Size and tube length for easy identification (Fig. 2)
- Transparent tube material for visual inspection of administration of nutrition and liquid
- Flexible, X-ray-opaque polyurethane olive with two openings at the distal tube tip for FR 8, FR 10, FR 12 (Fig. 3)
- Tube is numbered at 5 cm steps and marked at 1 cm steps
- Product package includes: Plaster for fixing the tube

	External ø	Internal ø	Length	gastric/ intestinal	Packaging	Art.-no.
Freka® Tubes, ENFit, Y-Connector						
Freka Tube FR 8, ENFit	2.8 mm	1.9 mm	80 cm	g	10 x 1	7981826
Freka Tube FR 8, ENFit	2.8 mm	1.9 mm	120 cm	g, i	10 x 1	7981827
Freka Tube FR 10, ENFit	3.4 mm	2.4 mm	120 cm	g, i	10 x 1	7981828
Freka Tube FR 12, ENFit	4.0 mm	2.9 mm	120 cm	g, i	10 x 1	7981829
Freka Tube FR 15, ENFit	4.5 mm	3.5 mm	120 cm	g, i	10 x 1	7981830

gastric /
intestinal

Freka Silicone



Freka Sil Tubes, ENFit are for gastric and intestinal feeding. Freka Sil Tubes, ENFit are made of silicone rubber.

- Enteral connector ENFit with linear connector (Fig. 1)
- X-ray-opaque tube material
- Tube markings every 10 cm
- Flexible, X-ray-opaque polyurethane olive with two openings at the distal tip for FR 8 and FR 13 tubes (Fig. 2)
- 5 lateral outlet holes (FR 16 tube)
- Product package includes: MCT oil as lubricant, adaptor for lubricant application and plaster for tube fixing

	External ø	Internal ø	Length	gastric/ intestinal	Packaging	Art.-no.
Freka® Tubes, ENFit, Silicone						
Freka Sil FR 8, ENFit	2.5 mm	1.6 mm	130 cm	g, i	10 x 1	7981835
Freka Sil FR 13, ENFit	4.3 mm	2.9 mm	120 cm	g, i	10 x 1	7981836
Freka Sil FR 16, ENFit	5.3 mm	3.5 mm	100 cm	g, i	10 x 1	7981837
Freka SIL-Sonde FR 16*	5.3 mm	3.5 mm	100 cm	g	10 x 1	7980951

*for decompression, with funnel connector

gastric /
intestinal